

## Enhancements in ZoomText 2020.2004.35 & Fusion 2020.2004.13 (April 2020)

The following changes are in the latest ZoomText 2020 and Fusion 2020 April releases. To download, select [ZoomText 2020](#) or [Fusion 2020](#). You must be logged on as an administrator to install this software.

- Resolved issues reported by the IRS related to an internal application running with ZoomText on Windows 10.
- Fixed font smoothing corruption that could happen when scrolling through Office 365 apps.
- Resolved an issue where dialog boxes were not being read correctly in Word 2016. (ZoomText only)
- Resolved an issue where the AppView and TextView features were not working in Word 2010. (ZoomText only)
- Resolved an issue where AppReader and Finder were not working with some apps on Windows 7 computers. (ZoomText only)
- Resolved an issue where Focus Enhancements, Cursor Enhancements, and mouse echo were not working with some apps on Windows 7 computers.
- Resolved an issue with ZoomText where mouse echo would return to and start reading from the top of a Word document while you are navigating by line through the document.
- Resolved an issue where ZoomText or Fusion would unexpectedly close. The problem appeared after installing recent Windows 10 updates.
- Resolved an issue where ZoomText would close after switching among apps (using **ALT+TAB**) and then changing magnification levels.
- Improved stability when running ZoomText with PowerPoint Office 365.
- Improved stability when running ZoomText with KDiff3 software.
- Addressed minor formatting issues in localized versions of the ZoomText 2020 user interface.

## Enhancements in JAWS 2020.2004.66 (April 2020)

The following is a list of improvements made between the March 2020 release and the April 2020 update.

### Announce Group Changes in Outlook 365 Message Lists

JAWS can now indicate group changes, such as dates or conversations, as you navigate through message lists. For example, if your Outlook view settings are set to show messages by date and you are navigating a folder containing messages covering multiple days, when you move to the first message in the new date group, you will hear announcements such as "today," "yesterday," "last week," or "two weeks ago." To turn on this feature, open Quick Settings (**INSERT+V**) while in Outlook and select the Announce Group Change in Lists check box.

### New Focus Braille Display Commands for Adjusting Speech Rate

You can now use the following controls on the Focus braille display to adjust the speech rate.

- Temporarily increase speech rate: **LEFT** or **RIGHT SHIFT+RIGHT SELECT CHORD**
- Temporarily decrease speech rate: **LEFT** or **RIGHT SHIFT+LEFT SELECT CHORD**
- Permanently increase speech rate: **RIGHT SELECT CHORD**
- Permanently decrease speech rate: **LEFT SELECT CHORD**

### Other Changes

- If Mouse Echo is enabled and you move over text with the mouse, JAWS or Fusion now only reads the current line as expected instead of reading everything on the screen.
- Book nine of the JAWS Basic Training, Surf's Up, Surfing the Internet with JAWS, an Introduction, has been updated to include links to go directly to the Surf's Up web pages. To get the latest book, select Training from the JAWS Help menu to open the Training Table of Contents in FSReader, select the link to open book nine, and then press **ENTER** when prompted to install the updated book.
- Added a new Quick Settings (**INSERT+V**) option called **List Nesting Level Announcement** which is available in Chrome and Edge Chromium. Use this option to toggle the announcement of nesting levels for list items in edit fields on web pages like Google Docs.
- In Microsoft Teams, when you press **CTRL+PERIOD** to display Teams hot keys, the Virtual Cursor is now active so you can review the keystrokes as if you were on a web page.
- Added updated braille drivers from Humanware containing performance improvements when using their line of displays with JAWS and Fusion.

## Download Links:

### **JAWS Downloads**

- JAWS Connected Link: [J2020.2004.66.400.exe](#)
- JAWS Offline Installers:
  - [J2020.2004.66.400-Offline-x64.exe](#)
  - [J2020.2004.66.400-Offline-x86.exe](#)

### **ZoomText Downloads**

- ZoomText Connected Link: [ZT2020.2004.35.400.exe](#)
- ZoomText Offline Installer: [ZT2020.2004.35.400-Offline.exe](#)

### **Fusion Downloads**

- Fusion Connected Link: [F2020.2004.13.400.exe](#)
- Fusion Offline Installer: [F2020.2004.13.400-Offline.exe](#)

## **Want to Learn More About Our Products and Services**

You can call our toll free number 888-211-6933 or email us at [info@nelowvision.com](mailto:info@nelowvision.com). Our two(2) state-of-the-art technology showrooms are a wonderful resource, showcasing a variety of leading electronic low vision magnifiers and blindness products. You can schedule an appointment with one of our Vision Specialists who can help you find the right technology, training and care solution.

If you cannot come to us, our staff can even drive out and bring our Vision Store to you, for a no-obligation product demonstration. We offer a comprehensive line of products including electronic magnification low vision aids, blindness products, and leading low vision products for people with glaucoma, macular degeneration and other forms of vision loss.