



# explorē 12

## Magnify your active lifestyle!

explorē 12 is an easy-to-use electronic magnifier with superb full high-definition image quality for improved near and distance viewing.

Its lightweight design makes it an ideal companion for on-the-go use, to maintain your active and mobile lifestyle.

As the largest portable magnifier in our **explorē** series, its 12-inch touchscreen is perfect for anyone desiring a larger interface.

When used with the foldable stand, **explorē 12** rises above the table, giving plenty of room to read and write.



Working height  
16.5cm / 6.5in

## Easier reading with live panning

The touchscreen's zoom in/out and panning functions allow you to read materials easily without having to move documents as much.



# Stay active and engaged in your world

With a superior magnification range, the tablet-sized **explorē 12** allows you to continue everyday activities like reading, writing, painting, viewing photos or maps, and much more.

magnification  
up to 30x



## Features

Polarized LED lighting to reduce glare from shiny surfaces

Distance camera

2-30x magnification with autofocus

Touchscreen gestures with adjustable brightness, auto shut-off, and preferred zoom

Image capture gallery

Available in 30 languages

17 color contrast options

2-year HumanWare warranty



## Specifications

12-inch, full high-definition screen with 2-integrated cameras

3.5-hour battery

Device weight: 2.5lbs / 1134g

Type-C USB charging and file transfer

High-resolution video output

## What else is in the box?

Folding, transportable reading stand (option)

Type-C USB and high-resolution video cables

Carry bag, neoprene protective sleeve, cleaning cloth and matte screen protector

User Manual and Getting Started (English/French)

**explorē 5, 8 and 12**  
Portable electronic magnifiers

**VictorReader Stratus**  
Talking book player

**Reveal 16i**  
Digital magnifier

## Low vision solutions for every need

When conventional glasses are insufficient and surgery or medical treatment is not appropriate, you should consider low vision aids.

Whether you experience macular degeneration, glaucoma, diabetic retinopathy, or cataracts, you can regain autonomy and re-engage with the world around you with a little help from HumanWare's solutions.

For more information on the **explorē 12**, please contact:

**USA** 1 800 722 3393 **CAN** 1 888 723 7273

or visit [humanware.com](http://humanware.com)

**Humanware**<sup>TM</sup>  
see things. differently.