

MagniLink iTAB



All-in-one solution for iPad

MagniLink iTAB is a portable and flexible solution that makes it possible to use an LVI video magnifier with your iPad. This combination becomes a versatile and useful visual aid in a compact format that is easy to take with you everywhere.

Compact solution for iPad

MagniLink iTAB is developed for and connects directly to Apple's iPad devices via the USB-C port. The ergonomic tiltable stand builds in comfort and an extra connector enables charging of the iPad while the system is in use. MagniLink iTAB can be ordered with a reading/distance camera with fantastic image quality. You can also connect an LVI video magnifier such as an MagniLink S*. Together with an iPad (not included), this becomes a versatile visual aid!

The flexible format is perfect for those who want to be able to use their system in different places, at home, in the workplace, or at school. A split function makes it possible to read simultaneously from a book or other printed text using the iPad's camera – and information on a whiteboard using iTAB's distance camera. It's also possible to split the screen with other applications on the iPad.

Software with limitless possibilities

The heart of the iTAB is our MagniLink Viewer software which is easily downloaded from the App Store. MagniLink Viewer offers different display modes that promote increased speed, endurance, and relaxation during reading. You can access digital and printed documents. The OCR function is always included, giving users the opportunity to have printed and digital text read aloud to stave off exhaustion. The MagniLink Viewer supports iPad's touch screen gestures, and short keys via a Bluetooth keyboard.

Capable combination

Together with the iPad apps, you get a very complete and capable system with multiple uses.

Alternative solutions for iPad

MagniWifi and MagniLink WifiCam are two solutions for those who want or need to be able to connect wirelessly to their iPad. Read more at www.lviglobal.se.





Model overview

Model*	Stand	Camera	Fits iPad-model**	Item nr
MagniLink iTAB SC 11	Yes	Yes	iPad Pro 11"	ML-iTAB-SC-11
MagniLink iTAB SC 12.9	Yes	Yes	lpad Pro 12.9"	ML-iTAB-SC-12.9
MagniLink iTAB S 11	Yes	No	lpad Pro 11"	ML-iTAB-S-11
MagniLink iTAB S 12.9	Yes	No	lpad Pro 12.9"	ML-iTAB-S-12.9

*iPad is not included. **There are specific system requirements for the iPad used with the MagniLink iTAB. Please see our website for updated system requirements (<https://lviamerica.com/catalog/products/magnilink-itable>). If MagniLink iTAB is purchased without camera, some MagniLink video magnifiers can be used instead, such as a MagniLink S, manufactured from 1st of November 2021.

Technical specifications

MagniLink iTAB

Weight 6.6 lb

Dimensions Folded: 14.4 x 11.8 x 3 In
Unfolded: 9 x 12.2 x 13.8 In



With camera



Without camera

External camera

Magnification 0.9-50x
(height 13", iPad Pro 12.9")

Frame rate 60 Hz

Length Camera and arm: 24.4 In

Weight 1.5 lb

MagniLink iTAB



LVI Low Vision International is recognized as one of the world's leading manufacturers of low vision aids with the mission to make daily life easier for people with visual impairments.

We are certified to ISO 14001 & 9001



Intertek

LVI America, Inc.

302 Saunders Road, Suite 200,
Riverwoods, IL 60015, USA
Phone: (888) 781-7811
E-mail: info@lviamerica.com

www.lviamerica.com



LVI Low Vision International
Verkstadsgatan 5
352 46 Växjö, SWEDEN
Tel: +46 (0)470 72 77 00
Fax: +46 (0)470 72 77 25
E-mail: info@lvi.se

WWW.LVI.SE